

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-8001

Weil-Blksly Rngr, Cer, str, 2.5mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Straight Cup Angle
- 200 mm Shaft Length

Weil-Blakesley Rongeur, Ceramic coated, straight, 2.5mm

http://buxtonbio.com/weilblksly_rngrcerstr2_5mm_886