

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-8201

Weil-BlkslyRngr, Cer, up30°2.5mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Up 30° Cup Angle
- 200 mm Shaft Length

Weil-Blakesley Rongeur, Ceramic coated, up 30°, 2.5mm

http://buxtonbio.com/weilblkslyrngrcerup302_5mm_888