

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-8261

W'l-Blaksly Rngr, Cer, up45°, 5mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Up 45° Cup Angle
- 200 mm Shaft Length

Weil-Blakesley Rongeur, Ceramic coated, up 45°, 5.0mm

http://buxtonbio.com/wlblaksly_rngrcerup455mm_890