

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-8522

Ferris-Smith Rngr, Cer, 180x2mm, dn



- Standard IVD
- Ceramic Black Finish
- Ferris-Smith Handle
- Down 30° Cup Angle

Ferris-Smith Rongeur, Ceramic coated, 180x2mm, down 30°

http://buxtonbio.com/frissmith_rngrcer180x2mmdn_913