

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone:973-560-4848
Fax: 973-560-4999

838-131-15

IVD Rngr,Cer,R-Hnd 130x1.5,str



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Right
- Straight Cup Angle
- Shaft Length 130

IVD Rongeur, Ceramic coated, right handled, 130mm shaft, 1.5mm cup straight

http://buxtonbio.com/ivd_rngrcerrhnd_130x1_5str_967