

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

838-162-15

IVD Rngr, Cer, R-Hnd 160x1.5, up30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Right
- Up 30° Cup Angle
- Shaft Length 160

IVD Rongeur, Ceramic coated, right handled, 160mm shaft, 1.5mm cup, up 30°

http://buxtonbio.com/ivd_rngrcerrhnd160x1_5up30_971