

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 839-112-20

IVD Rngr Cer L-hndl 110x2,up30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Left
- Up 30° Cup Angle
- Shaft Length 110

IVD Rongeur, Ceramic coated, left handled, 110mm shaft, 2mm cup, up 30°

[http://buxtonbio.com/ivd\\_rngr\\_cer\\_lhndl\\_110x2up30\\_959](http://buxtonbio.com/ivd_rngr_cer_lhndl_110x2up30_959)