

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

839-182-20

IVD Rngr, Cer, L-Hd, 180x2mm, up30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Left
- Up 30° Cup Angle
- Shaft Length 180

IVD Rongeur, Ceramic coated, left handled, 180mm shaft, 2mm cup, up 30°

http://buxtonbio.com/ivd_rngrcerlhd180x2mmup30_10097