

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

839-184-30

IVD Rngr, Cer, L-Hd, 180x3mm, dn30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Left
- Down 30° Cup Angle
- Shaft Length 180

IVD Rongeur, Ceramic coated, left handled, 180mm shaft, 3mm cup, down 30°

http://buxtonbio.com/ivd_rngrcerlhd180x3mmdn30_10098