

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

84-3563

Sm-Hndl Punch, 250x6mm, 90° dwn

Regular footplate
90° Downbiter



- Original Construction
- Down 90° Angled Foot Plate
- Lumbar; Regular Bite
- Small Handle
- 250 mm Shaft Length

Small-Handle Punch, 250x6mm, 90°
down-angled

http://buxtonbio.com/smhndl_punch_250x6mm90dwn_10081