

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 84-5041

Punch, 180x4mm, 90° up, TFP



- Original Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling None
- 180 mm Shaft Length

Sm-Handle Punch, 180x4mm, 90° up-angled, TFP

[http://buxtonbio.com/punch\\_180x4mm90uptfp\\_1064](http://buxtonbio.com/punch_180x4mm90uptfp_1064)