

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

84-5443

Punch, 230x4mm, 90° dwn, TFP

Thin footplate
90° Downbiter



- Original Construction
- Down 90° Angled Foot Plate
- Cervical; Regular Bite
- 230 mm Shaft Length

Sm-Hndl Punch, 230x4mm, 90° down-
angled, TFP

http://buxtonbio.com/punch_230x4mm90dwntfp_1109