

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-2062

LrgHnd Pnch, 180x6mm, 40° up, Ejec



- Larger 16 mm opening
- Ejector tip clears bone chips with each bite
- Up 40° Angled Foot Plate
- 180 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, 180x6mm, 40° up-angled, Ejec

http://buxtonbio.com/lrg_hnd_pnch180x6mm40upejec_1013