

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-2451

LrgHnd Pnch, 230x5mm, 90° up, Ejec



- Larger 16 mm opening
- Ejector tip clears bone chips with each bite
- Up 90° Angled Foot Plate
- 230 mm Shaft Length
- 5 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, 230x5mm, 90° up-angled, Ejec

http://buxtonbio.com/lrg_hnd_pnch230x5mm90upejec_1059