

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-2562

LrgHnd Pnch, 250x6mm, 40° up, Ejec



- Larger 16 mm opening
- Ejector tip clears bone chips with each bite
- Up 40° Angled Foot Plate
- 250 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, 250x6mm, 40° up-angled, Ejec

[http://buxtonbio.com/lrg\\_hnd\\_pnch250x6mm40upejec\\_1016](http://buxtonbio.com/lrg_hnd_pnch250x6mm40upejec_1016)