

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-3051

Punch, 180x5mm, 90° up, Ejec



- Original Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Small Handle
- 180 mm Shaft Length

Sm-Handle Punch, 180x5mm, 90° up-angled, Ejec

http://buxtonbio.com/punch_180x5mm90upejec_1048