

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-3532

Punch, 250x3mm, 40° up, Ejec



- Original Construction
- Up 40° Angled Foot Plate
- Lumbar; Regular Bite
- Chip Handling Ejector
- Small Handle
- 250 mm Shaft Length

Sm-Handle Punch, 250x3mm, 40° up-angled, Ejec

http://buxtonbio.com/punch_250x3mm40upejec_1007