

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-4451

Frs-Sm Pnch, 230x5mm, 90° up, Ejec



- Original Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Ferris-Smith Handle
- 230 mm Shaft Length

Frs-Smith Punch, 230x5mm, 90° up-angle, Ejec

http://buxtonbio.com/frssm_pnch230x5mm90upejec_1055