

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-4462

Frs-Sm Pnch, 230x6mm, 40° up, Ejec



- Original Construction
- Up 40° Angled Foot Plate
- Lumbar; Regular Bite
- Chip Handling Ejector
- Ferris-Smith Handle
- 230 mm Shaft Length

Frs-Smith Punch, 230x6mm, 40° up-angle, Ejec

[http://buxtonbio.com/frssm\\_pnch230x6mm40upejec\\_1011](http://buxtonbio.com/frssm_pnch230x6mm40upejec_1011)