

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-4941

Frs-Sm Pnch, 280x4mm, 90° up, Ejec



- Original Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Ferris-Smith Handle
- 280 mm Shaft Length

Frs-Smith Punch, 280x4mm, 90° up-angle, Ejec

http://buxtonbio.com/frssm_pnch280x4mm90upejec_1056