

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-7744

Bay-Kers'Up, RR, 40°, 180x4mm, TFP



Photo  
Pending

- Satin
- Original
- Offset Up
- Shaft Length 180mm
- Chip Clearing Ejector
- Foot Plate Style u40T

Bayonet Kerrison rongeur, 180 mm shaft, offset up, Ramrod ejector, 40° x4mm bite, thin footplate

[http://buxtonbio.com/baykersuprr40180x4mmtfp\\_1492](http://buxtonbio.com/baykersuprr40180x4mmtfp_1492)