

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-7754

Bay-Kers'Up, RR, 40°, 180x5mm, TFP



Photo
Pending

- Satin
- Original
- Offset Up
- Shaft Length 180mm
- Chip Clearing Ejector
- Foot Plate Style u40T

Bayonet Kerrison rongeur, 180 mm shaft, offset up, Ramrod ejector, 40° x5mm bite, thin footplate

http://buxtonbio.com/baykersuprr40180x5mmtfp_1492