

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-8132

Bay-Ker, R-hnd, RR, 110x3, 40°



- Satin
- Original
- Offset Right Handed
- Shaft Length 110mm
- Chip Set Ejector
- Up 40° Angled Foot Plate

Bayonet Kerrison rongeur, 110 mm shaft, right-handled, Ramrod ejector, 40° x3mm bite

http://buxtonbio.com/baykerrhndrr110x340_1530