

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-8222

Bay-Ker, R-hnd, RR, 150x2, 40°



Photo
Pending

- Satin
- Original
- Offset Right Handed
- Shaft Length 150mm
- Chip Set Ejector
- Foot Plate Angle Up 40°

Bayonet Kerrison rongeur, 150 mm shaft, right-handled, Ramrod ejector, 40° x2mm bite

http://buxtonbio.com/baykerrhndrr150x240_1538