

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-8244

Bay-Ker, R-hnd, RR, 150x4, 40°, TFP



- Satin
- Original
- Offset Right Handed
- Shaft Length 150mm
- Chip Set Ejector
- Up 40° Angled Thin Foot Plate

Bayonet Kerrison rongeur, 150 mm shaft, right-handled, Ramrod ejector, 40° x4mm bite, thin footplate

http://buxtonbio.com/baykerrhndrr150x440tfp_1539