

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-8341

Bay-Ker, R-hnd, RR, 180x4, 90°



- Satin
- Original
- Offset Right Handed
- Shaft Length 180mm
- Chip Set Ejector
- Up 90° Angled Foot Plate

Bayonet Kerrison rongeur, 180 mm shaft, right-handled, Ramrod ejector, 90° x4mm bite

http://buxtonbio.com/baykerrhndrr180x490_1546