

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-8361

Bay-Ker, R-hnd, RR, 180x6, 90°



Photo  
Pending

- Satin
- Original
- Offset Right Handed
- Shaft Length 180mm
- Chip Set Ejector
- Up 90° Angled Foot Plate

Bayonet Kerrison rongeur, 180 mm shaft, right-handled, Ramrod ejector, 90° x6mm bite

[http://buxtonbio.com/baykerrhndrr180x690\\_1546](http://buxtonbio.com/baykerrhndrr180x690_1546)