

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-8522

Bay-Ker, L-hnd, RR, 110x2, 40°



- Satin
- Original
- Offset Left Handed
- Shaft Length 110mm
- Chip Set Ejector
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, 110 mm shaft, left-handed, Ramrod ejector, 40° x2mm bite

http://buxtonbio.com/baykerlhndrr110x240_1502