

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-8622

Bay-Ker, L-hnd, RR, 150x2, 40°



Photo  
Pending

- Satin
- Original
- Offset Left Handed
- Shaft Length 150mm
- Chip Set Ejector
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, 150 mm shaft, left-handed, Ramrod ejector, 40° x2mm bite

[http://buxtonbio.com/baykerlhndrr150x240\\_1511](http://buxtonbio.com/baykerlhndrr150x240_1511)