

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-8724

Bay-Ker, L-hnd, RR, 180x2, 40°, TFP



- Satin
- Original
- Offset Left Handed
- Shaft Length 180mm
- Chip Set Ejector
- 40° Angled Thin Foot Plate

Bayonet Kerrison rongeur, 180 mm shaft, left-handed, Ramrod ejector, 40° x2mm bite, thin footplate

[http://buxtonbio.com/baykerlhndrr180x240tftp\\_1520](http://buxtonbio.com/baykerlhndrr180x240tftp_1520)