

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone:973-560-4848
Fax: 973-560-4999

845-5160

Foram.Rngr, lite curv,Set3,tfp

- Set of three in 2, 3 and 4 mm bites
- Foraminotomy punch with light "shallow" curve
- 200 mm shafts
- Original Construction
- Satin Finish
- Thin Footplates

Shallow Curve



Foraminotomy Rongeurs, light curve,
thin footplate, Set of 3, 2mm, 3mm,
4mm