

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-2043

LrgHnd Pnch, Cer, 180x4mm, 90°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 90° Angled Foot Plate
- 180 mm Shaft Length
- 4 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 180x4mm, 90°  
down-angle

[http://buxtonbio.com/lrgHnd\\_pnchcer180x4mm90dwn\\_1330](http://buxtonbio.com/lrgHnd_pnchcer180x4mm90dwn_1330)