

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-2231

LrgHnd Punch, Cer, 200x3mm, 90°up



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 90° Angled Foot Plate
- 200 mm Shaft Length
- 3 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 200x3mm, 90°  
up-angle

[http://buxtonbio.com/lrghnd\\_punchcer200x3mm90up\\_1264](http://buxtonbio.com/lrghnd_punchcer200x3mm90up_1264)