

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2561

Sm-Hdl Punch, Cer, 250x6mm, 90° up



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 90° Angled Foot Plate
- 250 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Sm-Hndl Punch, Ceramic, 250x6mm,
90° up-angle

http://buxtonbio.com/smhdl_punchcer250x6mm90up_1267