

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2963

LrgHnd Pnch, Cer, 300x6mm, 90°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 90° Angled Foot Plate
- 300 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 300x6mm, 90°
down-angle

http://buxtonbio.com/lrgHnd_pnchcer300x6mm90dwn_1334