

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-3412

### Sm-Hdl Punch, Cer, 230x1mm, 40°up

Regular footplate  
Ceramic Coated



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Lumbar; Regular Bite
- Chip Handling None
- Small Handle
- 230 mm Shaft Length

Sm-Hndl Punch, Ceramic, 230x1mm,  
40° up-angle

[http://buxtonbio.com/smhdl\\_punchcer230x1mm40up\\_1211](http://buxtonbio.com/smhdl_punchcer230x1mm40up_1211)