

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-4012

Frs-Sm Pnch, Cer, 180x2.5, 40° up



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Lumbar; Regular Bite
- Chip Handling None
- Ferris-Smith Handle
- 180 mm Shaft Length

Frs-Sm Punch, Cer, 180x2.5mm, 40°
up-angle

http://buxtonbio.com/frssm_pnchcer180x2_540up_1214