

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-4032

Frs-Sm Punch, Cer, 180x3mm, 40° up



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Lumbar; Regular Bite
- Chip Handling None
- Ferris-Smith Handle
- 180 mm Shaft Length

Frs-Sm Punch, Ceramic, 180x3mm,  
40° up-angle

[http://buxtonbio.com/frssm\\_punchcer180x3mm40up\\_1214](http://buxtonbio.com/frssm_punchcer180x3mm40up_1214)