

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-5062

Cer Kers Rngr, tfp, 50up, 180x6mm



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Cervical; Regular Bite
- Chip Handling None
- 180 mm Shaft Length

Kerrison Rongeur, Ceramic coated,  
thin footplate, 50°up, 180x6mm

[http://buxtonbio.com/cer\\_kers\\_rngrtfp50up180x6mm\\_1242](http://buxtonbio.com/cer_kers_rngrtfp50up180x6mm_1242)