

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-5064

Cer Kers Rngr, tfp, 50dn, 180x6mm



- Ceramic Black Finish
- 0
- Down 40° Angled Foot Plate
- Foot Plate Thickness Thin
- 180 mm Shaft Length

Kerrison Rongeur, Ceramic coated,
thin footplate, 50°down, 180x6mm

http://buxtonbio.com/cer_kers_rngrtfp50dn180x6mm_1313