

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-5422

Punch, Cer, 230x2mm, 40° up, TFP



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Cervical; Regular Bite
- Chip Handling None
- 230 mm Shaft Length

Sm-Hndl Punch, Cer, 230x2mm, 40°  
up-ang, TFP

[http://buxtonbio.com/punchcer230x2mm40uptfp\\_1244](http://buxtonbio.com/punchcer230x2mm40uptfp_1244)