

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-5432

Punch, Cer, 230x3mm, 40° up, TFP



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Cervical; Regular Bite
- Chip Handling None
- 230 mm Shaft Length

Sm-Hndl Punch, Cer, 230x3mm, 40°
up-ang, TFP

http://buxtonbio.com/punchcer230x3mm40uptfp_1244