

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-7460

X-ordinary IVD Rngr, Cer 30°up



- 300 mm shaft
- Ferris Smith ring handle
- Heavy-duty pattern
- Serrated jaws
- 8 mm bite
- Ceramic Coated
- Angled up 30°

Extraordinary IVD Rongeurs, Ceramic coated, angled 30°upwards

[http://buxtonbio.com/xordinary\\_ivd\\_rngr\\_cer\\_30up\\_920](http://buxtonbio.com/xordinary_ivd_rngr_cer_30up_920)