

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-7470

X-ordinary IVD Rngr, Cer 60°up



- 300 mm shaft
- Ferris Smith ring handle
- Heavy-duty pattern
- Serrated jaws
- 8 mm bite
- Ceramic Coated
- Angled up 60°

Extraordinary IVD Rongeurs, Ceramic coated, angled 60°upwards

http://buxtonbio.com/xordinary_ivd_rngr_cer_60up_920