

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-7621

Punch, Cer, U-Bay, 150x2mm, 90°



- Ceramic
- Original
- Offset Up
- Shaft Length 150mm
- Chip Set None
- 90° Angled Foot Plate

Sm-Hndl Punch, Cer, Up-Bayonet,
150x2mm, 90°

http://buxtonbio.com/punchcerubay_150x2mm90_1668