

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-7634

Pnch,Cer,U-Bay, 150x3mm,40°TFP

Photo  
Pending

- Ceramic
- Original
- Offset Up
- Shaft Length 150mm
- Chip Set None
- 40° Angled Thin Foot Plate

Sm-Hndl Punch, Cer,Up-Bay,  
150x3mm, 40°, TFP

[http://buxtonbio.com/pnchcerubay\\_150x3mm40tfp\\_1671](http://buxtonbio.com/pnchcerubay_150x3mm40tfp_1671)