

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-2062

LrgHnd Pnch, Cer, 180x6, 40°up, Ej



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Ejector Pin clears bone chips with each squeeze
- Up 40° Angled Foot Plate
- 180 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

LrgHndl Punch, Cer, 180x6mm, 40°
up-ang, Ejec

http://buxtonbio.com/lrgHnd_pnchcer180x640upej_1235