

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-2221

LrgHnd Pnch, Cer, 200x2, 90°up, Ej



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Ejector Pin clears bone chips with each squeeze
- Up 90° Angled Foot Plate
- 200 mm Shaft Length
- 2 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

LrgHndl Punch, Cer, 200x2mm, 90°  
up-ang, Ejec

[http://buxtonbio.com/lrgHnd\\_pnchcer200x290upej\\_1282](http://buxtonbio.com/lrgHnd_pnchcer200x290upej_1282)