

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-2241

LrgHnd Pnch, Cer, 200x4, 90°up, Ej



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Ejector Pin clears bone chips with each squeeze
- Up 90° Angled Foot Plate
- 200 mm Shaft Length
- 4 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

LrgHndI Punch, Cer, 200x4mm, 90°  
up-ang, Ejec

[http://buxtonbio.com/lrgHnd\\_pnchcer200x490upej\\_1282](http://buxtonbio.com/lrgHnd_pnchcer200x490upej_1282)