

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-2451

LrgHnd Pnch, Cer, 230x5, 90°up, Ej



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Ejector Pin clears bone chips with each squeeze
- Up 90° Angled Foot Plate
- 230 mm Shaft Length
- 5 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

LrgHndl Punch, Cer, 230x5mm, 90°
up-ang, Ejec

http://buxtonbio.com/lrghnd_pnchcer230x590upej_1283